

medel®

ELITE



Fully automatic
Clear
Comfortable

 **Soft inflate technology**



Extra-large backlit display



Irregular heartbeat indicator



Average of the last 3 readings



5-year warranty

Medel **ELITE** is a fully automatic arm blood pressure monitor; it is supplied with the innovative Soft Inflate Technology making measurement faster and more accurate. Inflation is soft and comfortable, allowing shorter measurement time and assuring the highest degree of accuracy.

Medel ELITE is supplied with T3 function that allows an average of the last 3 readings to be calculated simply and automatically. With a large backlit display for easy and clear measurement reading.

T3 function

it simply and automatically displays the average of the last 3 readings

Soft Inflate Technology

fully automatic, reliable and comfortable measuring system

ESH

tested in accordance with the ESH (European Society of Hypertension) protocol

ELITE

Medel ELITE is equipped with a universal cuff.
It is supplied with a sturdy case for travelling and transportation.



Technical specifications	ELITE
Technology	Oscillometric method
Accuracy	Pressure +/- 3 mmHg Heart rate +/- 5%
Measurement range	Pressure 0 to 300 mmHg; heart rate 40 to 199 beats/min.
Power supply	4 AA 1,5V batteries
Memory	60
Systolic/Diastolic and Pulse display	✓
Date and time	✓
Irregular heartbeat indicator (IHB)	✓
T3 Function	✓
LCD size	9,4 x 10,2 cm
Cuff Type	Universal 22-42 cm

COD. 92587 Medel ELITE

CE 0123

